



MEDIKAL-MUNERIS®
AL SERVICIO DE TU SALUD S.A DE C.V.

A common software platform for EEG, LTM, ICU, Sleep and research studies

What's new in NeuroWorks

- Windows® 10 compatibility
- Support for Quantum on-board memory and Nicolet Cortical Stimulator integrations
- Support for new cutting-edge Brain Monitor
- Frequency Analysis Tool
- Addition of Full HD (1080p) video

The Natus NeuroWorks platform simplifies the process of collecting, monitoring and managing data for routine EEG testing, long-term monitoring, ICU monitoring and Sleep studies. NeuroWorks systems are scalable to meet the needs of private practice clinics, hospitals and large teaching facilities.

Intuitive Clinical Experience

- Multiple Spike & Event analyzer options
- Fully synchronized HD video
- Monitoring options for home ambulatory studies
- Scalable remote access with Citrix provisioning
- Powerful qEEG trending

Streamlined Workflow

- Easy data searching with a variety of customizable parameters
- Automated settings synchronization across networked NeuroWorks stations
- Convenient import and export with EDF/EDF+

Time-Saving Administrator Functions

- Ease of setup using IP streaming cameras and amplifiers
- Powerful enterprise level SQL database
- Automated discovery of amplifiers and cameras using DHCP

natus®

NeuroWorks configures to your network and security needs

Robust Security Features

- HIPAA requirements via dedicated XLSecurity layer
- Role-based security access to features and studies
- Automated power recovery
- Microsoft® BitLocker™ drive encryption support

Advanced IT Capabilities

- Bi-directional HL7 integration with HIS/EMR
- Citrix® and VMWare® Ready partner – the most powerful infrastructure to facilitate collaboration



A cutting-edge, single solution for EEG, LTM, ICU and Sleep Studies

Advanced software for clinical excellence

- EEG/Sleep combo software package
- User customizable personal workspaces
- Software Development Kit for direct data access for research
- Data-Stream Wizard with two data streams for efficient clinical and research workflow*

Reliable and flexible hardware options

- Patient wearable compact and lightweight amplifiers
- Range of EEG amplifiers with 32- to 256- channels, plus differential and DC channel options
- Wireless amplifier options for long-term monitoring
- Range of sampling rates from 200 Hz to 16 kHz

A family of amplifiers to meet your needs



NEW! Natus Quantum



EMU40EX



NEW! Natus Brain Monitor



EEG32U



Trex HD



**DISTRIBUIDOR
EXCLUSIVO PARA
MÉXICO**

Neurology Service Programs o toda nuestra línea de suministros, póngase en contacto con su distribuidor

Número de teléfono para clientes de .:

33-38348777

ventas@medikal-muneris.com.mx

*Only available with Quantum amplifier

natus

Natus Medical Incorporated

natus.com