

# Comprehensive Solutions for Hospitals

EEG, LTM, ICU, EMG, IOM, Sleep & Supplies

- Empowering quality patient care
- Leading-edge equipment
- Advanced IT and connectivity solutions
- Unparalleled service and support



## Our Natus Brand Promise

Your partner in advancing patient care

As a global leader, Natus Neuro offers a **comprehensive portfolio** of **healthcare equipment and supplies** used by clinicians to **improve patient lives**.



## Natus Neuro Cybersecurity Program

Proactively addressing cybersecurity as an integral part of our design process

- Secure handling of confidential information generated or maintained by our hardware and software
- Continuous evaluation of our products against identified threats and vulnerabilities
- Monitoring of feedback channels to manage security events in the field

## Natus NeuroWorks® Solutions

A common software platform for EEG, LTM, ICU, Sleep and Research studies

### Leading-Edge Features


- Windows® 10 compatibility
- Intuitive user interface for simplified management of EEG studies
- Fully synchronized, easy to set up, Full HD (1080p) TCP/IP PTZ network cameras

### Advanced IT Capabilities & Remote Monitoring

- HL7 communication with HIS/EMR and Natus DataShare
- IT test environment system packages available
- Cerebrum Enterprise Solution – Automated software update rollout

CITRIX  
READY

vmware  
READY  
APPLICATION  
SOFTWARE



### Security & HIPAA Compliance

- XLSecurity Layer
- Role-based security package with drive encryption
- Automatic user activity logging

## A family of amplifiers to meet your needs:



Natus Quantum



Brain Monitor



EEG32U



Trex HD



EMU40EX



## Natus NeuroWorks LTM

An easy-to-use, state-of-the-art system for long-term EEG monitoring

- Automated data pruning
- Full HD (1080p) video
- A wide range of configuration options including:
  - Wall or ceiling mounted dual video systems
  - Compact, ergonomic and portable carts
  - Fixed installations
- Frequency analysis tool
- Multiple Spike & Event analyzer options
- Powerful trending options, including qEEG
- Stellate Gridview software option for assistance with functional brain mapping

CITRIX  
READY

vmware  
READY  
APPLICATION  
SOFTWARE

## Cortical Stimulator

Functional brain mapping for brain tumor and epilepsy surgery


- Use standalone or software controlled
- Displays selected and actual current delivered
- User configurable frequency, duration and current levels
- Ictal Disrupt feature for management of "after discharges"



## Natus Quantum® Amplifier

A Quantum Leap in amplifier technology in an ultra-small package

- Clinical EEG, long-term scalp and intracranial epilepsy monitoring and research
- 256 referential channels, up to 16 differential channels
- Excellent signal quality with up to 16 kHz sampling rate
- On-board memory with automated data backfill
- Integrated Digital Switch Matrix (DSM)



Software control integration in Natus NeuroWorks

# Routine EEG & ICU Monitoring Solutions

## Natus NeuroWorks EEG

Sophisticated EEG solution with  
an intuitive user interface

- On-board impedance check
- Integrated pulse oximeter
- Customized photic stimulator sequences
- Low noise signal acquisition
- True TCP/IP amplifiers and cameras with DHCP support plug anywhere on the network

**CITRIX**  
READY

**vmware**  
READY  
APPLICATION  
SOFTWARE



## Natus NeuroWorks ICU

Continuous ICU brain monitoring  
whenever and wherever you need it

- All-in-one touch-screen computer with intuitive UI designed for intensive workflow
- Flexible installation options - wall, column, boom or mobile cart
- Optional Neurowand for portable, low-profile, remote monitoring
- Powerful trending options, including qEEG
- Customized panel for nursing staff and advanced views for physicians

**CITRIX**  
READY

**vmware**  
READY  
APPLICATION  
SOFTWARE



Neurowand

Browse our full line of supplies  
and place orders:

[www.natusneurostore.com](http://www.natusneurostore.com)  
1-800-356-0007



## EEG and ICU Supplies

- Grass disposable and reusable cup electrodes
- Subdermal needle electrodes
- EEG Headcaps
- Creams, gels and pastes

## Trex HD and Trex HD Monitoring Systems

Reliable monitoring for physicians,  
freedom and flexibility for patients

### Trex HD System

Reliable, simple workflow for  
the patient home environment



- Up to 96 hours of HD video, wirelessly synchronized with EEG
- No cables or computer to set up at home – avoid costly patient home visits
- Lightest ambulatory amplifier in the industry
- Perfect for both EEG and full featured sleep studies

### Trex HD Monitoring System

Live monitoring – Remotely access and  
view data in real-time



## Protektor32 IOM System


Assurance, versatility, protection

- 16- or 32-channel multi-modality monitoring
- 16 high-level, 2 low-level programmable electrical stimulation outputs
- User-selectable electrode switching
- Detachable acquisition and stimulation pods



## IOM Supplies

- Disposable and reusable needle electrodes
- Corkscrew electrodes
- Disposable IOM probes





# EMG & Neuromuscular Ultrasound Solutions

## Nicolet® EDX EMG/NCS/EP/IOM System

Premium performance,  
maximum efficiency

- Multi-modality EMG/EP/NCS testing
- Advanced technology amplifiers provide superior data and reduced acquisition time
- Automatic markers, measurement tables and data storage
- Flexible data control with roll back, roll forward, averaging, auto store, signal enhancer
- Advanced testing options including Multi-MUP, CMAP, SFEMG and MUNIX
- 16-channel clinical EP and IOM option
- HL7 communication with HIS/EMR

**CITRIX®** **vmware®**



- Disposable adhesive electrodes
- TECA® Injectable needle electrodes
- Disposable and reusable concentric needle electrodes

## Natus UltraPro S100 EMG/NCS/EP System

Reliability, performance, portability  
– exceptional value

- Fast, intuitive clinical workflow
- Exceptional audio quality
- Auditory, visual and somatosensory EP options
- HL7 communication with HIS/EMR



## Natus InVisus Ultrasound System

Powerful, high-quality ultrasound in the palm of your hand

- Hand-held, HD tablet-based ultrasound system
- 5 – 18 MHz linear probe options
- Superior image quality for easy identification of structures and landmarks
- M-mode/B-mode, plus Color, Power, and Pulsed Wave Doppler
- Conduct needle insertions and injections with confidence using ultrasound guidance



## EMG Supplies

Browse our full line of supplies  
and place orders:

[www.natusneurostore.com](http://www.natusneurostore.com)  
1-800-356-0007

## Natus SleepWorks™ Software

Comprehensive, flexible and secure  
sleep diagnostic software

### Enhanced workflow for sleep studies

- Configurable Heads Up Display and observation notes
- Analyze and review ongoing studies, including video
- Customizable therapy types that simplify titration and reporting
- Robust security and site management features for enhanced productivity



Standard sleep testing  
with **Embla SDx**



Advanced sleep monitoring  
with **Embla NDx**

## Natus SleepWorks Amplifiers

Cutting-edge hardware that brings together  
long standing reliability and advanced technology

### Natus Base Unit Features

- Touchscreen display for real time interaction with ongoing studies
- Visual impedance checking and adding notes
- Perform bio-calibrations while in the patient room
- Universally compatible with several Natus EEG/PSG amplifier breakouts

### Breakout Box Features

- Small, rugged and lightweight
- Powered sensor inputs
- Integrated pressure sensor and pulse oximeter
- Up to 103 channels

## Embletta® MPR

The scalable sleep system that grows with your needs

- Integrated with Natus SleepWorks software
- Scalable from Level IV to Level I PSG
- Bio-Sensor Auto Start
- Multi-night recording for HST



Embletta MPR PG

ST+ Proxy

ST Proxy

## Sleep Supplies

- Respiratory effort sensors and belts
- Oximetry, airflow and body position sensors
- Disposable and reusable electrodes and snap on cables
- Creams, gels and pastes



# Service & Support



## Natus Neuro Training Academy

[www.natusneurologyacademy.com](http://www.natusneurologyacademy.com)

An online portal dedicated to educational programs and resources to support the use and optimization of Natus products.



Live and Recorded  
E-Seminars



In-House Clinical and  
IT/Biomedical Training



Comprehensive  
Clinical Library

## Neurodiagnostic Supplies

Convenient, complete, trusted

Natus Neuro supports the full spectrum of neurodiagnostic care, providing seamless solutions across all of our modalities.

See the Natus Supplies Catalog for a complete list of our products at [www.natus.com/supplies](http://www.natus.com/supplies)

US Customers visit our online store: [www.natusneurostore.com](http://www.natusneurostore.com)



**Natus Neuro**  
3150 Pleasant View Road  
Middleton, WI 53562 USA  
Tel: 1-800-356-0007  
1-608-829-8500  
Fax: 1-608-829-8709  
[www.natus.com](http://www.natus.com)

Para obtener más información sobre los productos de Natus, póngase en contacto con su distribuidor o representante de ventas local.

Número de teléfono para clientes de México.:  
**33-38348777**

[ventas@medikal-muneris.com.mx](mailto:ventas@medikal-muneris.com.mx)

## Natus Neuro Service & Support

A global service organization here for you

Natus Neuro Essential Support Elements

- Accessible and effective Technical Support
- Definitive technical documentation and knowledgeable installation teams
- Excellent replacement unit or spare part availability
- Extended warranty and service coverage programs



**DISTRIBUIDOR  
EXCLUSIVO PARA  
MÉXICO**

